

PLXNB3 Polyclonal Antibody

Catalog No: #29703



Package Size: #29703-1 50ul #29703-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PLXNB3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human PLXNB3 (NP_001156729.1).
Other Names	PLXNB3; PLEXB3; PLEXR; PLXN6; plexin-B3
Accession No.	Swiss-Prot#:Q9ULL4NCBI Gene ID:5365
Uniprot	Q9ULL4
GeneID	5365;
Calculated MW	170kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse brain, using PLXNB3 at 1:1000 dilution.

Background

The protein encoded by this gene is a member of the plexin family. It functions as a receptor for semaphorin 5A, and plays a role in axon guidance, invasive growth and cell migration. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only