

SNRPD3 Polyclonal Antibody

Catalog No: #29706



Package Size: #29706-1 50ul #29706-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SNRPD3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SNRPD3 (NP_001265585).
Other Names	SNRPD3; SMD3; Sm-D3; small nuclear ribonucleoprotein Sm D3
Accession No.	Swiss-Prot#:P62318NCBI Gene ID:6634
Uniprot	P62318
GeneID	6634;
Calculated MW	14kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using SNRPD3 at 1:1000 dilution.

Background

This gene encodes a core component of the spliceosome, which is a nuclear ribonucleoprotein complex that functions in pre-mRNA splicing. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only