

MLF2 Polyclonal Antibody

Catalog No: #29710



Package Size: #29710-1 50ul #29710-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MLF2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MLF2 (NP_005430.1).
Other Names	MLF2; NTN4; myeloid leukemia factor 2
Accession No.	Swiss-Prot#:Q15773NCBI Gene ID:8079
Uniprot	Q15773
GeneID	8079;
Calculated MW	28kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:100 - 1:200


Images




Immunohistochemistry of paraffin-embedded rat liver using MLF2 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human esophageal using MLF2 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse spleen using MLF2 at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using MLF2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only