

TMEM59 Polyclonal Antibody

Catalog No: #29719



Package Size: #29719-1 50ul #29719-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TMEM59 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human TMEM59 (NP_004863.2).
Other Names	TMEM59; C1orf8; DCF1; HSPC001; PRO195; UNQ169; transmembrane protein 59
Accession No.	Swiss-Prot#:Q9BXS4NCBI Gene ID:9528
Uniprot	Q9BXS4
GeneID	9528;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

This gene encodes a protein shown to regulate autophagy in response to bacterial infection. This protein may also regulate the retention of amyloid precursor protein (APP) in the Golgi apparatus through its control of APP glycosylation. Overexpression of this protein has been found to promote apoptosis in a glioma cell line. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only