

INTS7 Polyclonal Antibody

Catalog No: #29729



Package Size: #29729-1 50ul #29729-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	INTS7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human INTS7 (NP_001186740.1).
Other Names	INTS7; C1orf73; INT7; integrator complex subunit 7
Accession No.	Swiss-Prot#:Q9NVH2NCBI Gene ID:25896
Uniprot	Q9NVH2
GeneID	25896;
Calculated MW	107kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using INTS7 at 1:1000 dilution.

Background

This gene encodes a subunit of the integrator complex. The integrator complex associates with the C-terminal domain of RNA polymerase II and mediates 3'-end processing of the small nuclear RNAs U1 and U2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Note: This product is for in vitro research use only