

UNC50 Polyclonal Antibody

Catalog No: #29730



Package Size: #29730-1 50ul #29730-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	UNC50 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human UNC50 (NP_054763.2).
Other Names	UNC50; GMH1; HSD23; PDLs22; UNCL; URP; protein unc-50 homolog
Accession No.	Swiss-Prot#:Q53HI1NCBI Gene ID:25972
Uniprot	Q53HI1
GeneID	25972;
Calculated MW	30kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using UNC50 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only