

ZNF777 Polyclonal Antibody

Catalog No: #29731



Package Size: #29731-1 50ul #29731-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF777 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZNF777 (NP_056509.2).
Other Names	ZNF777; zinc finger protein 777
Accession No.	Swiss-Prot#:Q9ULD5NCBI Gene ID:27153
Uniprot	Q9ULD5
GeneID	27153;
Calculated MW	85kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF777 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only