

TMEM176B Polyclonal Antibody

Catalog No: #29733



Package Size: #29733-1 50ul #29733-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TMEM176B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TMEM176B (NP_054739.3).
Other Names	TMEM176B; LR8; MS4B2; transmembrane protein 176B
Accession No.	Swiss-Prot#:Q3YBM2NCBI Gene ID:28959
Uniprot	Q3YBM2
GeneID	28959;
Calculated MW	29kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TMEM176B at 1:1000 dilution.

Background

Note: This product is for in vitro research use only