

PDZRN4 Polyclonal Antibody

Catalog No: #29734



Package Size: #29734-1 50ul #29734-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PDZRN4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PDZRN4 (NP_037509.3).
Other Names	PDZRN4; LNX4; SAMCAP3L; PDZ domain containing ring finger 4
Accession No.	Swiss-Prot#:Q6ZMN7NCBI Gene ID:29951
Uniprot	Q6ZMN7
GeneID	29951;
Calculated MW	117kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PDZRN4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only