

HN1 Polyclonal Antibody

Catalog No: #29738



Package Size: #29738-1 50ul #29738-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	HN1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human HN1 (NP_057269.1).
Other Names	JPT1; ARM2; HN1; HN1A; Jupiter microtubule associated homolog 1
Accession No.	Swiss-Prot#:Q9UK76NCBI Gene ID:51155
Uniprot	Q9UK76
GeneID	51155;
Calculated MW	16kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of HT-29 cells, using HN1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only