

ARL15 Polyclonal Antibody

Catalog No: #29739

Package Size: #29739-1 50ul #29739-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARL15 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ARL15 (NP_061960.1).
Other Names	ARL15; ARFRP2; ADP ribosylation factor like GTPase 15
Accession No.	Swiss-Prot#:Q9NXU5NCBI Gene ID:54622
Uniprot	Q9NXU5
GeneID	54622;
Calculated MW	23kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse pancreas, using ARL15 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only