

STAP2 Polyclonal Antibody

Catalog No: #29741



Package Size: #29741-1 50ul #29741-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	STAP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human STAP2 (NP_001013863).
Other Names	STAP2; BKS; signal-transducing adaptor protein 2
Accession No.	Swiss-Prot#:Q9UGK3NCBI Gene ID:55620
Uniprot	Q9UGK3
GeneID	55620;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse heart, using STAP2 at 1:1000 dilution.



Western blot analysis of extracts of HT-29 cells, using STAP2 at 1:1000 dilution.

Background

This gene encodes the substrate of breast tumor kinase, an Src-type non-receptor tyrosine kinase. The encoded protein possesses domains and several tyrosine phosphorylation sites characteristic of adaptor proteins that mediate the interactions linking proteins involved in signal transduction pathways. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only