

ZNF701 Polyclonal Antibody

Catalog No: #29742



Package Size: #29742-1 50ul #29742-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF701 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ZNF701 (NP_060730.2).
Other Names	ZNF701; zinc finger protein 701
Accession No.	Swiss-Prot#:Q9NV72NCBI Gene ID:55762
Uniprot	Q9NV72
GeneID	55762;
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF701 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only