

CHPT1 Polyclonal Antibody

Catalog No: #29743



Package Size: #29743-1 50ul #29743-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CHPT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human CHPT1 (NP_064629.2).
Other Names	CHPT1; CPT; CPT1; cholinephosphotransferase 1
Accession No.	Swiss-Prot#:Q8WUD6NCBI Gene ID:56994
Uniprot	Q8WUD6
GeneID	56994;
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of DU145 cells, using CHPT1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only