

EBF4 Polyclonal Antibody

Catalog No: #29744

Package Size: #29744-1 50ul #29744-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EBF4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human EBF4 (NP_001103984.1).
Other Names	EBF4; COE4; O/E-4; early B-cell factor 4
Accession No.	Swiss-Prot#:Q9BQW3NCBI Gene ID:57593
Uniprot	Q9BQW3
GeneID	57593;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

EBF4 belongs to the conserved Olf/EBF family of helix-loop-helix transcription factors, members of which play important roles in neural development and B-cell maturation.

Note: This product is for in vitro research use only