

CHID1 Polyclonal Antibody

Catalog No: #29746

Package Size: #29746-1 50ul #29746-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CHID1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human CHID1 (NP_076436.3).
Other Names	CHID1; GL008; SI-CLP; SICLP; chitinase domain containing 1
Accession No.	Swiss-Prot#:Q9BWS9NCBI Gene ID:66005
Uniprot	Q9BWS9
GeneID	66005;
Calculated MW	55kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse pancreas, using CHID1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only