

TMC5 Polyclonal Antibody

Catalog No: #29749



Package Size: #29749-1 50ul #29749-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TMC5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human TMC5 (NP_079056.2).
Other Names	TMC5; transmembrane channel like 5
Accession No.	Swiss-Prot#:Q6UXY8NCBI Gene ID:79838
Uniprot	Q6UXY8
GeneID	79838;
Calculated MW	114kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TMC5 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only