

WBSCR27 Polyclonal Antibody

Catalog No: #29761

Package Size: #29761-1 50ul #29761-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	WBSCR27 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human WBSCR27 (NP_689772.2).
Other Names	METTL27;WBSCR27
Accession No.	Uniprot:Q8N6F8GeneID:155368
Uniprot	Q8N6F8
GeneID	155368
Calculated MW	26kDa
SDS-PAGE MW	26kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded Human colon using WBSCR27 Rabbit pAb.



Western blot analysis of extracts of various cell lines, using WBSCR27 antibody.



Immunohistochemistry of paraffin-embedded Rat testis using WBSCR27 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Mouse brain using WBSCR27 Rabbit pAb.

Background

This gene encodes a protein belonging to ubiE/COQ5 methyltransferase family. The gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.22-q11.23.

Note: This product is for in vitro research use only