

CCDC110 Polyclonal Antibody

Catalog No: #29763



Package Size: #29763-1 50ul #29763-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CCDC110 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CCDC110 (NP_689988.1).
Other Names	CCDC110; CT52; KM-HN-1; KMHN1; coiled-coil domain containing 110
Accession No.	Swiss-Prot#:Q8TBZ0NCBI Gene ID:256309
Uniprot	Q8TBZ0
GeneID	256309;
Calculated MW	97kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse heart, using CCDC110 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only