

YIPF7 Polyclonal Antibody

Catalog No: #29764

Package Size: #29764-1 50ul #29764-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	YIPF7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human YIPF7 (NP_872398.2).
Other Names	YIPF7; FinGER9; protein YIPF7
Accession No.	Swiss-Prot#:Q8N8F6NCBI Gene ID:285525
Uniprot	Q8N8F6
GeneID	285525;
Calculated MW	31kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using YIPF7 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only