

ZNF677 Polyclonal Antibody

Catalog No: #29765



Package Size: #29765-1 50ul #29765-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF677 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human ZNF677 (NP_872415.1).
Other Names	ZNF677; zinc finger protein 677
Accession No.	Swiss-Prot#:Q86XU0NCBI Gene ID:342926
Uniprot	Q86XU0
GeneID	342926;
Calculated MW	68kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF677 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only