

GTPBP8 Polyclonal Antibody

Catalog No: #29768



Package Size: #29768-1 50ul #29768-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GTPBP8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GTPBP8 (NP_054889.2).
Other Names	GTPBP8;HSPC135
Accession No.	Uniprot:Q8N3Z3GeneID:29083
Uniprot	Q8N3Z3
GeneID	29083
Calculated MW	Refer to figures
SDS-PAGE MW	32kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using GTPBP8 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Human liver using GTPBP8 Rabbit pAb.



Immunofluorescence analysis of L929 cells using GTPBP8 antibody.



Immunohistochemistry of paraffin-embedded Rat kidney using GTPBP8 Rabbit pAb.

Note: This product is for in vitro research use only