

GET4 Polyclonal Antibody

Catalog No: #29769



Package Size: #29769-1 50ul #29769-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GET4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GET4 (NP_057033.2).
Other Names	GET4; C7orf20; CEE; CGI-20; TRC35; golgi to ER traffic protein 4
Accession No.	Swiss-Prot#:Q7L5D6NCBI Gene ID:51608
Uniprot	Q7L5D6
GeneID	51608;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using GET4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only