

MAGOHB Polyclonal Antibody

Catalog No: #29770

Package Size: #29770-1 50ul #29770-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MAGOHB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human MAGOHB (NP_060518.1).
Other Names	MAGOHB;MGN2;mago;magoh
Accession No.	Uniprot:Q96A72GeneID:55110
Uniprot	Q96A72
GeneID	55110
Calculated MW	Refer to figures
SDS-PAGE MW	17kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using MAGOHB Rabbit pAb.

Background

Required for pre-mRNA splicing as component of the spliceosome. Plays a redundant role with MAGOH in the exon junction complex and in the nonsense-mediated decay (NMD) pathway.

Note: This product is for in vitro research use only