

RSPH10B Polyclonal Antibody

Catalog No: #29771



Package Size: #29771-1 50ul #29771-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RSPH10B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human RSPH10B (NP_775836.3).
Other Names	RSPH10B; radial spoke head 10 homolog B
Accession No.	Swiss-Prot#:P0C881NCBI Gene ID:222967
Uniprot	P0C881
GeneID	222967;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only