

ZNF620 Polyclonal Antibody

Catalog No: #29772

Package Size: #29772-1 50ul #29772-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF620 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ZNF620 (NP_001243096.1).
Other Names	ZNF620; zinc finger protein 620
Accession No.	Swiss-Prot#:Q6ZNG0NCBI Gene ID:253639
Uniprot	Q6ZNG0
GeneID	253639;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse thymus, using ZNF620 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only