

SDHAF2 Polyclonal Antibody

Catalog No: #29776

Package Size: #29776-1 50ul #29776-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SDHAF2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human SDHAF2 (NP_060311.1).
Other Names	SDHAF2; C11orf79; PGL2; SDH5; succinate dehydrogenase complex assembly factor 2
Accession No.	Swiss-Prot#:Q9NX18NCBI Gene ID:54949
Uniprot	Q9NX18
GeneID	54949;
Calculated MW	15kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SDHAF2 at 1:1000 dilution.

Background

This gene encodes a mitochondrial protein needed for the flavination of a succinate dehydrogenase complex subunit required for activity of the complex. Mutations in this gene are associated with paraganglioma.

Note: This product is for in vitro research use only