

GSPT1 Polyclonal Antibody

Catalog No: #29779



Package Size: #29779-1 50ul #29779-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GSPT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GSPT1 (NP_001123478.1).
Other Names	GSPT1; 551G9.2; ETF3A; GST1; eRF3a; G1 to S phase transition 1
Accession No.	Swiss-Prot#:P15170NCBI Gene ID:2935
Uniprot	P15170
GeneID	2935;
Calculated MW	80kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using GSPT1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only