

VASN Polyclonal Antibody

Catalog No: #29780

Package Size: #29780-1 50ul #29780-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VASN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human VASN (NP_612449.2).
Other Names	VASN; SLITL2; vasorin
Accession No.	Swiss-Prot#:Q6EMK4NCBI Gene ID:114990
Uniprot	Q6EMK4
GeneID	114990;
Calculated MW	71kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using VASN at 1:1000 dilution.

Background

Note: This product is for in vitro research use only