

DPPA5 Polyclonal Antibody

Catalog No: #29794



Package Size: #29794-1 50ul #29794-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DPPA5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DPPA5 (NP_001020461.1).
Other Names	DPPA5; ESG1; developmental pluripotency associated 5
Accession No.	Swiss-Prot#:A6NC42NCBI Gene ID:340168
Uniprot	A6NC42
GeneID	340168;
Calculated MW	13kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DPPA5 at 1:1000 dilution.

Background

This gene encodes a protein that may function in the control of cell pluripotency and early embryogenesis. Expression of this gene is a specific marker for pluripotent stem cells. Pseudogenes of this gene are located on the short arm of chromosome 10 and the long arm of chromosomes 14 and 19.

Note: This product is for in vitro research use only