

## FOXQ1 Polyclonal Antibody

Catalog No: #29796

Package Size: #29796-1 50ul #29796-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FOXQ1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	A synthetic peptide of human FOXQ1
Other Names	FOXQ1; HFH1; forkhead box Q1
Accession No.	Swiss-Prot#:Q9C009NCBI Gene ID:94234
Uniprot	Q9C009
GeneID	94234;
Calculated MW	41kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of rat kidney, using FOXQ1 at 1:1000 dilution.

## Background

FOXQ1 is a member of the FOX gene family, which is characterized by a conserved 110-amino acid DNA-binding motif called the forkhead or winged helix domain. FOX genes are involved in embryonic development, cell cycle regulation, tissue-specific gene expression, cell signaling, and tumorigenesis.

---

Note: This product is for in vitro research use only