

MYCL Polyclonal Antibody

Catalog No: #29801



Package Size: #29801-1 50ul #29801-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MYCL Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MYCL (NP_005367.2).
Other Names	MYCL; L-Myc; LMYC; MYCL1; bHLHe38; protein L-Myc
Accession No.	Swiss-Prot#:P12524NCBI Gene ID:4610
Uniprot	P12524
GeneID	4610;
Calculated MW	50kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using MYCL at 1:1000 dilution.



Western blot analysis of extracts of A-549 cells, using MYCL at 1:1000 dilution.

Background

Note: This product is for in vitro research use only