

PAG1 Polyclonal Antibody

Catalog No: #29802

Package Size: #29802-1 50ul #29802-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PAG1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	A synthetic peptide of human PAG1 (NP_060910.3).
Other Names	PAG1; CBP; PAG; phosphoprotein membrane anchor with glycosphingolipid microdomains 1
Accession No.	Swiss-Prot#:Q9NWQ8NCBI Gene ID:55824
Uniprot	Q9NWQ8
GeneID	55824;
Calculated MW	46kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PAG1 at 1:1000 dilution.

Background

The protein encoded by this gene is a type III transmembrane adaptor protein that binds to the tyrosine kinase csk protein. It is thought to be involved in the regulation of T cell activation.

Note: This product is for in vitro research use only