

DUSP15 Polyclonal Antibody

Catalog No: #29806



Package Size: #29806-1 50ul #29806-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DUSP15 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DUSP15 (NP_542178.2).
Other Names	DUSP15; C20orf57; VHY; dual specificity phosphatase 15
Accession No.	Swiss-Prot#:Q9H1R2NCBI Gene ID:128853
Uniprot	Q9H1R2
GeneID	128853;
Calculated MW	32kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DUSP15 at 1:1000 dilution.

Background

The protein encoded by this gene has both protein-tyrosine phosphatase activity and serine/threonine-specific phosphatase activity, and therefore is known as a dual specificity phosphatase. This protein may function in the differentiation of oligodendrocytes. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only