

MYRF Polyclonal Antibody

Catalog No: #29817



Package Size: #29817-1 50ul #29817-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | MYRF Polyclonal Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Affinity purification |
| Applications | WB |
| Species Reactivity | Human,Mouse |
| Immunogen Description | A synthetic peptide of human MYRF (NP_001120864.1). |
| Other Names | MYRF; C11orf9; MRF; Ndt80; pqn-47; myelin regulatory factor |
| Accession No. | Swiss-Prot#:Q9Y2G1NCBI Gene ID:745 |
| Uniprot | Q9Y2G1 |
| GeneID | 745; |
| Calculated MW | Refer to figures |
| Formulation | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage | Store at -20°C |

Application Details

WB□1:500 - 1:2000

Background

This gene encodes a transcription factor that is required for central nervous system myelination and may regulate oligodendrocyte differentiation. It is thought to act by increasing the expression of genes that effect myelin production but may also directly promote myelin gene expression. Loss of a similar gene in mouse models results in severe demyelination. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only