

COL4A5 Polyclonal Antibody

Catalog No: #29819



Package Size: #29819-1 50ul #29819-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	COL4A5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human COL4A5 (NP_203699.1).
Other Names	COL4A5; ASLN; ATS; CA54; collagen alpha-5(IV) chain
Accession No.	Swiss-Prot#:P29400NCBI Gene ID:1287
Uniprot	P29400
GeneID	1287;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IHC□1:50 - 1:200

Background

This gene encodes one of the six subunits of type IV collagen, the major structural component of basement membranes. Mutations in this gene are associated with X-linked Alport syndrome, also known as hereditary nephritis. Like the other members of the type IV collagen gene family, this gene is organized in a head-to-head conformation with another type IV collagen gene so that each gene pair shares a common promoter. Alternatively spliced transcript variants have been identified for this gene.

Note: This product is for in vitro research use only