

Slc31a2 Polyclonal Antibody

Catalog No: #29820

Package Size: #29820-1 50ul #29820-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Slc31a2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	A synthetic peptide of mouse Slc31a2 (NP_001277447.1).
Other Names	SLC31A2; COPT2; CTR2; hCTR2; AI604396; solute carrier family 31 member 2
Accession No.	Swiss-Prot#:O15432NCBI Gene ID:20530
Uniprot	O15432
GeneID	1318;
Calculated MW	16kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse pancreas, using Slc31a2 at 1:500 dilution.

Background

Note: This product is for in vitro research use only