

GRINA Polyclonal Antibody

Catalog No: #29826



Package Size: #29826-1 50ul #29826-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	GRINA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human GRINA (NP_000828.1).
Other Names	GRINA; HNRGW; LFG1; NMDARA1; TMBIM3; protein lifeguard 1
Accession No.	Swiss-Prot#:Q7Z429NCBI Gene ID:2907
Uniprot	Q7Z429
GeneID	2907;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

IHC 1:100 - 1:200

Images



Western blot analysis of extracts of Mouse pancreas, using GRINA antibody at 1:1000 dilution.Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.Detection: ECL Basic Kit.Exposure time: 1s.

Note: This product is for in vitro research use only