

FOXA3 Polyclonal Antibody

Catalog No: #29827



Package Size: #29827-1 50ul #29827-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FOXA3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	A synthetic peptide of human FOXA3 (NP_004488.2).
Other Names	FOXA3; FKHH3; HNF3G; TCF3G; forkhead box A3
Accession No.	Swiss-Prot#:P55318NCBI Gene ID:3171
Uniprot	P55318
GeneID	3171;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. The crystal structure of a similar protein in rat has been resolved.

Note: This product is for in vitro research use only