

NUBP1 Polyclonal Antibody

Catalog No: #29830

Package Size: #29830-1 50ul #29830-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUBP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NUBP1 (NP_002475.2).
Other Names	NUBP1; NBP; NBP1; NBP35; nucleotide binding protein 1
Accession No.	Swiss-Prot#:P53384NCBI Gene ID:4682
Uniprot	P53384
GeneID	4682;
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded rat heart using NUBP1 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human thyroid cancer using NUBP1 at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human colon carcinoma using NUBP1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse liver using NUBP1 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using NUBP1 at dilution of 1:100 (40x lens).



Western blot analysis of extracts of HepG2 cells, using NUBP1 at 1:1000 dilution.

Background

NUBP1 is a member of the NUBP/MRP subfamily of ATP-binding proteins.

Note: This product is for in vitro research use only