

NUBP1 Polyclonal Antibody

Catalog No: #29831



Package Size: #29831-1 50ul #29831-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NUBP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human NUBP1 (NP_002475.2).
Other Names	NUBP1; NBP; NBP1; NBP35; nucleotide binding protein 1
Accession No.	Swiss-Prot#:P53384NCBI Gene ID:4682
Uniprot	P53384
GeneID	4682;
Calculated MW	35kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NUBP1 at 1:1000 dilution.

Background

NUBP1 is a member of the NUBP/MRP subfamily of ATP-binding proteins.

Note: This product is for in vitro research use only