

ZNF182 Polyclonal Antibody

Catalog No: #29842



Package Size: #29842-1 50ul #29842-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF182 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF182 (NP_001171570.1).
Other Names	ZNF182; HHZ150; KOX14; ZNF21; Zfp182; zinc finger protein 182
Accession No.	Swiss-Prot#:P17025NCBI Gene ID:7569
Uniprot	P17025
GeneID	7569;
Calculated MW	74kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF182 at 1:1000 dilution.

Background

Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This gene encodes a zinc finger protein, and belongs to the krueppel C2H2-type zinc-finger protein family. It may be involved in transcriptional regulation. Multiple alternatively spliced transcript variants encoding different isoforms have been identified.

Note: This product is for in vitro research use only