

PLAA Polyclonal Antibody

Catalog No: #29846



Package Size: #29846-1 50ul #29846-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PLAA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PLAA (NP_001026859.1).
Other Names	PLAA; DOA1; PLA2P; PLAP; phospholipase A2 activating protein
Accession No.	Swiss-Prot#:Q9Y263NCBI Gene ID:9373
Uniprot	Q9Y263
GeneID	9373;
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded rat ovary using PLAA at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human placenta using PLAA at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using PLAA at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only