

CAND2 Polyclonal Antibody

Catalog No: #29855

Package Size: #29855-1 50ul #29855-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CAND2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human CAND2 (NP_001155971.1).
Other Names	CAND2;TIP120B;Tp120b
Accession No.	Uniprot:O75155GeneID:23066
Uniprot	O75155
GeneID	23066
Calculated MW	130kDa
SDS-PAGE MW	130kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded Mouse spleen using CAND2 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Human colon carcinoma using CAND2 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Rat testis using CAND2 Rabbit pAb.



Western blot analysis of extracts of various cell lines, using CAND2 antibody.

Background

Probable assembly factor of SCF (SKP1-CUL1-F-box protein E3 ubiquitin ligase complexes that promotes the exchange of the substrate-recognition F-box subunit in SCF complexes, thereby playing a key role in the cellular repertoire of SCF complexes.

Note: This product is for in vitro research use only