

MON2 Polyclonal Antibody

Catalog No: #29856



Package Size: #29856-1 50ul #29856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MON2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human MON2 (NP_001265399.1).
Other Names	MON2; protein MON2 homolog
Accession No.	Swiss-Prot#:Q7Z3U7NCBI Gene ID:23041
Uniprot	Q7Z3U7
GeneID	23041;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Background

Note: This product is for in vitro research use only