

FSCN2 Polyclonal Antibody

Catalog No: #29857

Package Size: #29857-1 50ul #29857-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FSCN2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human FSCN2 (NP_036550.1).
Other Names	FSCN2; RFSN; RP30; fascin-2
Accession No.	Swiss-Prot#:O14926NCBI Gene ID:25794
Uniprot	O14926
GeneID	25794;
Calculated MW	55kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using FSCN2 at 1:1000 dilution.

Background

This gene encodes a member of the fascin protein family. Fascins crosslink actin into filamentous bundles within dynamic cell extensions. This family member is proposed to play a role in photoreceptor disk morphogenesis. A mutation in this gene results in one form of autosomal dominant retinitis pigmentosa and macular degeneration. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only