

DCAF12 Polyclonal Antibody

Catalog No: #29858



Package Size: #29858-1 50ul #29858-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DCAF12 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human DCAF12 (NP_056212.1).
Other Names	DCAF12; CT102; KIAA1892; TCC52; WDR40A; DDB1 and CUL4 associated factor 12
Accession No.	Swiss-Prot#:Q5T6F0NCBI Gene ID:25853
Uniprot	Q5T6F0
GeneID	25853;
Calculated MW	47kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DCAF12 at 1:1000 dilution.

Background

This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates.

Note: This product is for in vitro research use only