

MYRIP Polyclonal Antibody

Catalog No: #29859



Package Size: #29859-1 50ul #29859-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MYRIP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human MYRIP (NP_001271354.1).
Other Names	MYRIP; SLAC2-C; SLAC2C; rab effector MyRIP
Accession No.	Swiss-Prot#:Q8NFW9NCBI Gene ID:25924
Uniprot	Q8NFW9
GeneID	25924;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only