

COMMD5 Polyclonal Antibody

Catalog No: #29862



Package Size: #29862-1 50ul #29862-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	COMMD5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human COMMD5 (NP_054785.2).
Other Names	COMMD5; HCARG; HT002; COMM domain containing 5
Accession No.	Swiss-Prot#:Q9GZQ3NCBI Gene ID:28991
Uniprot	Q9GZQ3
GeneID	28991;
Calculated MW	25kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using COMMD5 at 1:1000 dilution.



Immunofluorescence analysis of C6 cells using COMMD5 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only