

# POP5 Polyclonal Antibody

Catalog No: #29867



Package Size: #29867-1 50ul #29867-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                            |
|-----------------------|----------------------------------------------------------------------------|
| Product Name          | POP5 Polyclonal Antibody                                                   |
| Host Species          | Rabbit                                                                     |
| Clonality             | Polyclonal                                                                 |
| Isotype               | IgG                                                                        |
| Purification          | Affinity purification                                                      |
| Applications          | WB                                                                         |
| Species Reactivity    | Human,Rat                                                                  |
| Immunogen Description | A synthetic peptide of human POP5 (NP_057002.2).                           |
| Other Names           | POP5; HSPC004; RPP2; RPP20; hPop5; ribonuclease P/MRP protein subunit POP5 |
| Accession No.         | Swiss-Prot#:Q969H6NCBI Gene ID:51367                                       |
| Uniprot               | Q969H6                                                                     |
| GeneID                | 51367;                                                                     |
| Calculated MW         | 19kDa                                                                      |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.          |
| Storage               | Store at -20°C                                                             |

## Application Details

WB 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using POP5 at 1:1000 dilution.

## Background

---

Note: This product is for in vitro research use only